

Code No. 18647

Anti-Human**14-3-3 γ Protein Rabbit IgG Affinity purify**

Volume : 100 μ g

-
- Introduction** : The 14-3-3 proteins are a family of conserved regulatory molecules expressed in all eukaryotic cells. A striking feature of the 14-3-3 proteins is their ability to bind a multitude of functionally diverse signaling proteins, including kinases, phosphatases, and transmembrane receptors. This plethora of interacting proteins allows 14-3-3 to play important roles in a wide range of vital regulatory processes, such as mitogenic signal transduction, apoptotic cell death, and cell cycle control.
- Antigen** : Synthetic peptide for C-terminal of Human 14-3-3 γ (QQDDDGEGNN)
- Purification** : Affinity Purified with antigen peptide
- Form** : Lyophilized product from 1% BSA in PBS containing 0.05%NaN₃
- How to use** : 1.0 mL deionized water will be added to the product
- Dilution** : PBS (pH7.4) containing 1% BSA
- Stability** : Lyophilized product, 5 years at 2 - 8 °C
: Solution, 2 years at -20 °C
- Application** : This antibody can be used for immunohistochemistry with formalin fixed paraffin embedded tissues after microwave treatment (for 10 mn, 10 mM citrate buffer, pH 6.0). The optimal concentration is about 5 μ g/mL, however, the concentration should be optimized by each laboratory.
: This antibody can be used for western blotting in concentration of 0.2 - 0.5 μ g/mL.
- Specificity** : Human 14-3-3 gamma specific
Non-crossreact to Human 14-3-3 β , ϵ , ζ , η , τ
- Reference** : 1. Hermeking H, Lengauer C, Polyak K, He TC, Zhang L, Thiagalingam S, Kinzler KW, Vogelstein B. 14-3-3 sigma is a p53-regulated inhibitor of G2/M progression. Mol Cell. 1997 Dec; 1(1): 3-11.
2. Chan TA, Hermeking H, Lengauer C, Kinzler KW, Vogelstein B. 14-3-3 Sigma is required to prevent mitotic catastrophe after DNA damage. Nature. 1999 Oct 7; 401(6753): 616-20.
3. Nakajima T, Shimooka H, Weixa P, Segawa A, Motegi A, Jian Z, Masuda N, Ide M, Sano T, Oyama T, Tsukagoshi H, Hamanaka K, Maeda M. Immunohistochemical demonstration of 14-3-3 sigma protein in normal human tissues and lung cancers, and the preponderance of its strong expression in epithelial cells of squamous cell lineage. Pathol Int. 2003 Jun; 53(6): 353-60.

For research use only, not for use in diagnostic procedures.

Distributed by:



Immuno-Biological Laboratories, Inc. Toll-Free: 888-523-1246

8201 Central Ave NE, Suite P

Minneapolis, MN 55432

Email: info@IBL-America.com

Web: www.IBL-America.com