



Code No. 18651

**Anti-Human  
Thioredoxin (C) Rabbit IgG Affinity Purify**

Volume : 100 µg

---

<b>Antigen</b>	: Synthetic peptide in portion of C-terminus of Human Thioredoxin
<b>Source</b>	: Serum
<b>Subclass</b>	: IgG
<b>Purification</b>	: Purified with antigen peptide
<b>Form</b>	: Lyophilized product from 1% BSA in PBS containing 0.05% NaN <sub>3</sub>
<b>How to use</b>	: 1.0 ml distilled water will be added to the product
<b>Dilution</b>	: PBS (pH7.4) containing 1% BSA
<b>Stability</b>	: Lyophilized product, 5 years at 2 – 8 °C : Solution, 2 years at –20 °C
<b>Application</b>	: This antibody can be used for western blotting in concentration of 2 ~ 5 µg/ml.

---

*For research use only, not for use in diagnostic procedures.*

Distributed by:



Immuno-Biological Laboratories, Inc. Toll-Free: 888-523-1246  
8201 Central Ave NE, Suite P Email: [info@IBL-America.com](mailto:info@IBL-America.com)  
Minneapolis, MN 55432 Web: [www.IBL-America.com](http://www.IBL-America.com)