

Product Information



Recombinant human TL1A/TNFSF15 protein

Catalog Number: IBATGP1049

PRODUCT INFORMATION

Expression system

E.coli

Domain

72-251aa

UniProt No.

O95150

NCBI Accession No.

NP_005109

Alternative Names

Tumor necrosis factor ligand superfamily member 15, TL1, TL1A, VEGI, VEGI192A

PRODUCT SPECIFICATION

Molecular Weight

22.7 kDa (201aa) confirmed by MALDI-TOF

Concentration

0.5mg/ml (determined by Bradford assay)

Formulation

Liquid in. 20mM Tris-HCl buffer (pH 8.5) containing 0.15M NaCl, 40% glycerol

Purity

> 90% by SDS-PAGE

Endotoxin level

< 1 EU per 1ug of protein (determined by LAL method)

Biological Activity

Measured by its binding ability in a functional ELISA with Human DR3/TNFRSF25 (CAT# ATGP4017)

Tag

His-Tag

Application

SDS-PAGE, Bioactivity

Storage Condition

Can be stored at +2C to +8C for 1 week. For long term storage, aliquot and store at -20C to -80C. Avoid repeated freezing and thawing cycles.

BACKGROUND

Description

Tumor necrosis factor ligand superfamily member 15, also known TNFSF15, belongs to the TNF superfamily.

For research use only. This product is not intended or approved for human, diagnostics or veterinary use.



Manufactured for:

Immuno-Biological Laboratories, Inc. (IBL-America)

8201 Central Ave. NE, Suite P, Minneapolis, Minnesota 55432, USA

Phone: (888) 523-1246

Fax: (763) 780-2988

Email: info@ibl-america.com

Web: www.ibl-america.com

Product Information



Recombinant human TL1A/TNFSF15 protein

Catalog Number: IBATGP1049

TNFSF15 is predominantly expressed in endothelial cells and its expression is inducible by TNFalpha and IL1alpha. TNFSF15 is an angiogenesis inhibitor of the TNF family and functions in part by directly inhibiting endothelial cell proliferation. TNFSF15 may act as an autocrine factor to induce apoptosis in endothelial cells via activation of multiple signaling pathways, including stress protein kinases and certain caspases. Recombinant human TNFSF15 protein, fused to His-tag at N-terminus, was expressed in E. coli and purified by using conventional chromatography techniques.

Amino acid Sequence

<MGSSHHHHHH SGLVPRGSH M>LKGQEFAPS HQQVYAPLRA DGDKPRALHT VVRQTPTQHF KNQFPALHWE
HELGLAFTKN RMNYTNKFLI IPESGDYFIY SQVTFRGMTS ECSEIRQAGR PNKPDSTVW ITKVTDSTPE PTQLLMGTKS
VCEVGSNWFQ PIYLGAMFSL QEGDKLMVNV SDISLVDYTK EDKTFFGAFL L


General References

Migone T S., et al. (2002) Immunity. 16:479.

Zhai Y., et al. (1999) FASEB J.13:181-189.

DATA


SDS-PAGE



3ug by SDS-PAGE under reducing condition and visualized by coomassie blue stain.

Biological Activity

Human TL1A/TNFSF15 is coated at 10 ug/ml (100 ul/well) can bind Human DR3/TNFRSF25 (CAT# ATGP4017) in a Functional ELISA assay.



For research use only. This product is not intended or approved for human, diagnostics or veterinary use.



Manufactured for:
Immuno-Biological Laboratories, Inc. (IBL-America)
8201 Central Ave. NE, Suite P, Minneapolis, Minnesota 55432, USA
Phone: (888) 523-1246 Fax: (763) 780-2988
Email: info@ibl-america.com Web: www.ibl-america.com