

Product information



INHA, 233-366aa


Human, His-tagged, Recombinant, *E.coli*

Cat. No. IBATGP1898

Full name: inhibin alpha chain

NCBI Accession No.: NP_002182

Description: INHA, also known as Inhibin alpha, belongs to the TGF-beta family. Inhibins and activins inhibit and activate, respectively, the secretion of follitropin by the pituitary gland. The inhibin alpha subunit joins either the beta A or beta B subunit to form a pituitary FSH secretion inhibitor. Inhibin has been shown to regulate gonadal stromal cell proliferation negatively and to have tumour-suppressor activity. Inhibins/activins are involved in regulating a number of diverse functions such as hypothalamic and pituitary hormone secretion, gonadal hormone secretion, germ cell development and maturation, erythroid differentiation, insulin secretion, nerve cell survival, embryonic axial development or bone growth, depending on their subunit composition. Inhibins appear to oppose the functions of activins. Recombinant human INHA protein, fused to His-tag at N-terminus, was expressed in *E.coli*.

<p>Form: Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.4M Urea, 10% glycerol</p> <p>Molecular Weight: 17.0kDa (155aa)</p> <p>Purity: > 85% by SDS - PAGE</p> <p>Concentration: 0.5 mg/ml (determined by Bradford assay)</p>	 <p>15% SDS-PAGE (3ug)</p>
<p>Sequences of amino acids:</p> <p><u>MGSSHHHHH</u> <u>SSGLVPRGSH</u> <u>MSTPLMSWPW</u> SPSALRLLQR PPEEPAAHAN CHRVALNLSF QELGWERWIV YPPSIFHYC HGGCGLHIPP NLSLPVPGAP PTPAQPYSLI PGAQPCCAAL PGTMRPLHVR TSDGGYSFK YETVPNLLTQ HCACI</p>	

General references:

Knight, P.G. (1996) *Front. Neuroendocrinol.* 17: 476-509.

Mather, J.P., *et al.* (1997). *Proc. Soc. Exp. Biol. Med.* 215: 209-222.

Storage: Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C.
Avoid repeated freezing and thawing cycles.

For research use only. This product is not intended or approved for human, diagnostics or veterinary use.



Manufactured for:

Immuno-Biological Laboratories, Inc. (IBL-America)
8201 Central Ave. NE, Suite P, Minneapolis, Minnesota 55432, USA
Phone: (888) 523-1246 Fax.: (763) 780-2988
Email: info@ibl-america.com Web: www.ibl-america.com