

Product information

RWDD4, 1-188aa

Human, His-tagged, Recombinant, *E.coli*

Cat. No. IBATGP2713

Full name: RWD domain-containing protein 4

NCBI Accession No.: NP_689895

Synonyms: RWD domain containing 4, FAM28A, RWDD4A


Description: RWDD4 is RWD domain-containing protein 4. Recombinant human RWDD4 protein, fused to His-tag at N-terminus, was expressed in *E.coli* and purified by using conventional chromatography techniques.

Form: Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 20% glycerol, 1mM DTT

Molecular Weight: 23.6 kDa (211aa) confirmed by MALDI-TOF

Purity: > 90% by SDS - PAGE

Concentration: 0.5 mg/ml (determined by Bradford assay)



Sequences of amino acids:

MGSSHHHHH SSGLVPRGSH MGSMSANEDQ EMELEALRSI YEGDESFREL SPVSFQYRIG ENGDPKAFLI EISWTETYPQ TPPILSMNAF
FNNTISSAVK QSI LAKLQEA VEANLGTAMT YTLFEYAKDN KEQFMENHNP INSATSISNI ISLETPTAP SSKKKDKKEQ LSKAQKRKLA
DKTDHKGELP RGWNWVDVVK HLSKTGSKDD E

General references:

Ota T., Suzuki Y., *et al.* (2004) *Nat. Genet.* 36:40-45

Storage: Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.

For research use only. This product is not intended or approved for human, diagnostics or veterinary use.



Manufactured for:

Immuno-Biological Laboratories, Inc. (IBL-America)
8201 Central Ave. NE, Suite P, Minneapolis, Minnesota 55432, USA
Phone: (888) 523-1246 Fax.: (763) 780-2988
Email: info@ibl-america.com Web: www.ibl-america.com