

Recombinant human AMPD2 protein

Catalog Number: IBATGP2875



PRODUCT INFORMATION

Expression system

E.coli

Domain

236-879aa

UniProt No.

Q01433

NCBI Accession No.

NP_004028.3

Alternative Names

AMP deaminase 2, adenosine monophosphate deaminase 2

PRODUCT SPECIFICATION

Molecular Weight

77 kDa (667aa)

Concentration

0.25mg/ml (determined by Bradford assay)

Formulation

Liquid in. Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Purity

> 95% by SDS-PAGE

Tag

His-Tag

Application

SDS-PAGE

Storage Condition

Can be stored at +2C to +8C for 1 week. For long term storage, aliquot and store at -20C to -80C. Avoid repeated freezing and thawing cycles.

BACKGROUND

Description

AMPD2 is important in purine metabolism by converting AMP to IMP. This protein, which acts as a homotetramer, is one of three AMP deaminases found in mammals. Several transcript variants encoding different isoforms have been found for this gene. Recombinant human AMPD2 protein, fused to His-tag at N-terminus, was expressed in E. coli and purified by using conventional chromatography techniques.

Amino acid Sequence

<MGSSHHHHHH SSGLVPRGSH MGS>DLLDAAK SVVRALFIRE KYMALSLQSF CPTTRRYLQQ LAEKPLETRT

For research use only. This product is not intended or approved for human, diagnostics or veterinary use.

Manufactured for:

Immuno-Biological Laboratories, Inc. (IBL-America)
8201 Central Ave. NE, Suite P, Minneapolis, Minnesota 55432, USA
Phone: (888) 523-1246 Fax.: (763) 780-2988
Email: info@ibl-america.com Web: www.ibl-america.com

Recombinant human AMPD2 protein

Catalog Number: IBATGP2875



YEQGPDPVPS ADAPVHPPAL EQHPYEHCEP STMPGDLGLG LRMVRGVVHV YTRREPDEHC SEVELPYPDL QEFVADVNVL
MALIINGPIK SFCYRRLQYL SSKFQMHVLL NEMKELAAQK KVPHRDFYNI RKVDTHIHAS SCMNQKHLLR FIKRAMKRHL
EEIVHVEQGR EQTLREVFES MNLTAYDLSV DTLDVHADRN TFHRFDKFNA KYNPIGESVL REIFIKTNR VSGKYFAHII
KEVMSDLEES KYQNAELRLS IYGRSRDEWD KLARWAVMHR VHSPNVRWL VQVPRLFDVYR TKGQLANFQE MLENIFLPLF
EATVHPASHP ELHLFLEHVD GFDSVDDESK PENHVFNL ES PLPEAWVEED NPPYAYYLYY TFANMAMLNH LRRQRGFHTF
VLRPHCGEAG PIHHLVSAFM LAENISHGLL LRKAPVLQYL YYLAQIGIAM SPLSNNSLFL SYHRNPLPEY LSRGLMVSL S
TDDPLQFHFT KEPLMEEYSI ATQVWKLSSC DMCELARNSV LMSGFSHKVK SHWLGPNYTK EGPEGNDIRR TNVPDIRVGY
RYETLCQELA LITQAVQSEM LETIPEEAGI TMSPGPQ


General References

Burkard T.R. et al. (2011) BMC Syst. Biol. 5:17-17.

Rigbolt K.T. et al. (2011) Sci. Signal. 4:RS3-RS3.

DATA

SDS-PAGE



3ug by SDS-PAGE under reducing condition and visualized by coomassie blue stain

For research use only. This product is not intended or approved for human, diagnostics or veterinary use.

Manufactured for:

Immuno-Biological Laboratories, Inc. (IBL-America)

8201 Central Ave. NE, Suite P, Minneapolis, Minnesota 55432, USA

Phone: (888) 523-1246

Fax.: (763) 780-2988

Email: info@ibl-america.com

Web: www.ibl-america.com